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## NOTE:

Power frequency magnetic field (50 Hz)	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical residential or commercial environment.
--	-------	-------	--

### RADIOFREQUENCY IMMUNITY

Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic environment - guidance
Conducted radio frequency	3 Veff from 150 kHz to	3 Veff from 150 kHz to	Portable and mobile RF communications equipment should be used no closer to any part of the <i>FR II POWER</i> device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance  $d = 1.2 \cdot \sqrt{P}$ from 150 kHz to 80 MHz  $d = 1.2 \cdot \sqrt{P}$ from 80 MHz to 800 MHz  $d = 2.3 \cdot \sqrt{P}$ from 800 MHz to 2.5 GHz  where <b>P</b> is the maximum output power rating of the transmitter in watts <b>W</b> according to the transmitter manufacturer and <b>d</b> is the separation distance in meters
Radiated radio frequency	3 Veff from 80 MHz to 2.5GHz	3 Veff from 80 MHz to 2.5GHz	

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Electrostatic discharge	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 6$ kV contact $\pm 8$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst	$\pm 2$ kV phase(s)-ground	$\pm 2$ kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Surge	$\pm 2$ kV phase(s)-ground $\pm 1$ kV	$\pm 2$ kV common mode $\pm 1$ kV differential	Mains power quality should be that of a typical residential or commercial environment.
Voltage dips, short interruptions and voltage variations on power supply input lines	< 5% $U_t$ (95% dip) for 0.5 cycles 40% $U_t$ (60% dip) for 5 cycles 70% $U_t$ (30% dip) for 25 cycles < 5% $U_t$ (95% dip) for 5 seconds	< 5% $U_t$ (95% dip) for 0.5 cycles 40% $U_t$ (60% dip) for 5 cycles 70% $U_t$ (30% dip) for 25 cycles < 5% $U_t$ (95% dip) for 5 seconds	Mains power quality should be that of a typical residential or commercial environment. If the user of the <i>FRII POWER</i> device requires continued operation during power mains interruptions, it is recommended that the <i>FRII POWER</i> device be powered from an uninterruptible power supply.

<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF Emissions CISPR 11	Group 1	The <i>FRII POWER</i> device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The <i>FRII POWER</i> device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies. Buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	The <i>FRII POWER</i> device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

<b>ELECTROMAGNETIC IMMUNITY</b>			
<b>Immunity Test</b>	<b>IEC 60601-1-2 Test Level</b>	<b>Compliance Level</b>	<b>Electromagnetic environment - guidance</b>

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

Problem description	Detail	Intervention
After having pressed button A "SPEED", the RED LED is steady-on.		System not operating. Replace the device with a new one.
The RED LED is steady-on.		Flat battery. Replace the central unit with a charged central unit.
After having pressed button A "SPEED" to increase the speed, the blue LEDs do not increase.		System not operating. The red LED should be lit. Replace the device with a new one.
After having pressed button B "CONTINUOUS ROTATION" or C "ALTERNATING ROTATION" in order to select the rotation mode, no modification of the rotation mode is seen.		System not operating. The red LED should be lit. Replace the device with a new one.


## 11. INFORMATION PROVIDED BY THE MANUFACTURER REGARDING THE ELECTROMAGNETIC COMPATIBILITY OF THE DEVICE

The information tables shown below have been compiled as required in point 6 of reference regulation EN 60601-1-2.

### ELECTROMAGNETIC EMISSIONS

The *FRII POWER* device is intended for use in the electromagnetic environment specified below. The customer or user of the *FRII POWER* device should assure that it is used in such an environment.

For more information on recycling, contact the Minister of the Environment or local authorities.

 **WARNING:** Do not dispose of this product as household waste. Arrange for the re-use or recycling of this product as set forth in Directive 2002/96/EC of the European Parliament and European Union Council on the disposal of waste electrical and electronic equipment (WEEE). The directive does not apply to contaminated products.

## 9. GUARANTEE

The following components of the device are guaranteed, only if used, stored and maintained correctly, for two years from the date of delivery:

- INTERNAL CENTRAL UNIT;
- CHARGING BASE.

The connection cables and all the single-use parts are not covered by the guarantee.

After this period, it is recommended that the device undergo periodic testing by the manufacturer every 2 years.

## 10. TROUBLESHOOTING

This section describes the most common problems that may arise during use of the *FRII Power* system.

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<b>Charging Base Power Supply</b>	110/230Vac ±15%
<b>Power Supply Frequency</b>	50/60 Hz
<b>Battery life</b>	1 Hour (full charge)
<b>Charging time</b>	2 Hours
<b>Insulation Class</b>	II
<b>Temperature of Use</b>	+5 - +40°C
<b>Humidity</b>	30% - 65%
<b>Controls</b>	3 buttons (A, B, C) and regulator
<b>Central Unit Dimensions</b>	170 x 23 x 36 mm
<b>Charging Base Dimensions</b>	260 x 140 x 230 mm
<b>Tool Body Dimensions</b>	365 x 55 x 38 mm
<b>Tool Body Weight including the central unit</b>	1070g

## 7. SERVICE LIFE OF THE INSTRUMENT

The *FRII Power* system is built to be long-lasting if the conditions of use and maintenance specified in this User Manual are met. Nevertheless, the service life of the complete device is considered to be the following:

- 2 Years for the INTERNAL CENTRAL UNIT;
- 3 Years for the CHARGING BASE;
- Single use for the TOOL BODY, which is provided in sterile and disposable form.

## 8. DISPOSAL

The *FRII Power* system contains electronic parts that must be treated as electronic waste. Dispose of the device and accessories in compliance with local regulations.

Abide by national recycling laws and the hospital's current recycling procedures in order to ensure proper disposal of the *FRII Power* device.

Any malfunction must be reported immediately to the manufacturer who will provide for repair or special maintenance. Only the manufacturer is authorized to carry out repairs or maintenance.

Where requested and with a valid Department of Clinical Engineering, according to our instructions and agreed estimates, any repairs or maintenance work may also be delegated to the aforesaid Department along with all related liabilities, Any tampering with the device, or any of its components, may invalidate correct operation of the device and void the existing guarantees.

## 5. PRESERVATION, STORAGE AND TRANSPORT

The *FRII Power* device must be used, stored and shipped in the following conditions:

<b>Storage temperature:</b>	<b>from 10°C to +40°C Maximum</b>
<b>degree of relative humidity:</b>	<b>65%</b>
<b>Atmospheric pressure:</b>	<b>from 700 hPa to 1060 hPa</b>

The CHARGING BASE must be stored in a dry location, protected against direct sunlight and any accidental water splashes.

Carefully avoid any type of impact for each component of the device.

Store in a secure area as specified in the "Warnings" section of this manual.

## 6. TECHNICAL INFORMATION

The *FRII Power* device is a battery-powered system designed in compliance with medical standards that guarantee patient safety.

Table 1 provides the technical specifications of the *FRII Power* system.

<b>Model</b>	<i>FRII Power</i>
<b>Classification</b>	Battery-powered system
<b>Casing Degree of Protection</b>	IP54
<b>Container</b>	Plastic
<b>Tool body Handpiece Power Supply</b>	Rechargeable 1.2 V battery

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

## WARNING

- The red-colored LED, when lit, indicates a device malfunction.
- In Continuous Rotation mode, do not exceed 12 minutes of continued use. For each 12 minutes of continuous use, switch the device off for 5 minutes.

## 3. CLEANING AND DISINFECTION

The CHARGING BASE can be cleaned taking great care not to damage it or allow any substances to enter the openings. Either "detergent" or "disinfectant" type liquids can be used, although care should be taken not to rub so hard as to damage the surfaces, especially those where the commands and symbols are located. It is recommended that soft fabric cloths be used which leave no residue or, better yet, disposable TNT cloths. The amount of liquid used must not be too much or such as to penetrate inside the container, since this could damage or interfere with the electrical parts inside. If this should occur, the CHARGING BASE has to be sent to a Clinical Engineering Department, if available at the user's facilities or structure, or to the manufacturer for the appropriate overhaul.

The CENTRAL UNIT must be cleaned using detergents and/or disinfectants normally used in hospitals and healthcare facilities, taking care not to rub so hard as to damage the surfaces. Soft brushes or sponges should be used, avoiding abrasive material completely in order to prevent damage.

The TOOL BODY is supplied in a sterile package, thus it must not be resterilized or reused after its first use.

## 4. MAINTENANCE

The CENTRAL UNIT does not require any routine maintenance; only periodic replacement of the rechargeable battery is required. The battery's expected life is 500 charge cycles. Depending on the frequency of use of the device, this period may range from 6 months to a year.

- To operate the device in “**ALTERNATING**” mode, press button **C**.

- To increase or decrease the speed of the device, after having executed the step described in point 9, press button **A**

- Each time the button is pressed, the graduated LED scale located on the body of the device increases as indicated by a series of lit green LEDs. The more LEDs that are lit, the higher the set speed.  
To reduce the speed, continue to press the button until the LEDs start over again from 0

- To stop the device, press the button corresponding to the set operating mode.

- If the “**CONTINUOUS ROTATION**” function is set, press button **B**,

- If the “**ALTERNATING**” function is set, press button **C**.

- To adjust the suction, use the curved slider located on the handpiece.

- Pushing the slider forward increases the suction flow, while pulling it back reduces the suction flow.

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


## 2.6 SUCTION

The device allows suction to be performed during the working cycle. In order for this function to be possible, the rear conical attachment on the TOOL BODY device must be connected to the centralized suction system of the operating room. The tube must be sterile.

The suction slider located on the upper part of the grip can be used to adjust the flow rate of suctioned water. Note that the flow rate is highly conditioned by the following factors, at the same rpm and power of the pump being used:

- the flow rate is reduced by increasing the rpm of the instrument;
- the flow rate is reduced if the channel is partially obstructed by residue of excised tissue;
- the flow rate is reduced if the suction tube is connected incorrectly, twisted or damaged.

## 2.7 INSTRUCTIONS FOR USE OF THE DEVICE

- Connect the connector assembled at the back of the device to the centralized suction system.

- Press button **A** "SPEED"  for 5 seconds to place the device in Stand-By.

- To operate the device in "CONTINUOUS ROTATION" mode, press button **B**.

- After use, the TOOL BODY device must be disposed of as special hospital waste according to current regulations and hospital instructions.

## 2.5 OPERATING PARAMETERS

With 1 green LED lit in continuous rotation, the speed is 500 RPM with each additional LED lit, the speed increases by 500 RPM/LED with a maximum of 6000 RPM.

With 1 green LED lit in alternating rotation, the speed is 500 RPM. With each additional LED lit, the speed increases by 500 RPM/LED with a maximum of 3000 RPM.


The red-colored LED, when lit, indicates a device malfunction.

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

UNIT from the TOOL BODY device, just open the END CAP by pressing lightly on the black rubber stripes on the side of the cap itself.

- Rotate the device by 45/60 degrees to allow the Central Unit to fall out by gravity.

- The CENTRAL UNIT can be cleaned using disinfectants normally used in hospitals and healthcare facilities, although care should be taken not to rub too hard as that could ruin the surfaces. Soft brushes or sponges should be used, avoiding abrasive material completely in order to prevent damage. The TOOL BODY is supplied in a sterile package, thus it must not be resterilized or reused after its first use.

- Non-sterile personnel can recover the Central Unit and replace it in the appropriate seat on the charging base.

- Personnel in non-sterile gowns will insert the CENTRAL UNIT inside the TOOL BODY with a firm maneuver, paying utmost attention to ensure that sterile and non-sterile parts do not come into contact with one another.

- Remove the protective cone inserted by the personnel in sterile gowns. Now close the end cap of the FRII POWER device.

- The device is now ready for operation.

- The blades / burs of the FRII POWER device can be rotated in 4 different positions. Unscrew the ring nut counterclockwise, pulling slightly until the blade/bur comes out of its seat. Position the blade in the desired position and firmly re-tighten the ring nut. Now the FRII POWER device is ready to be used.

- At the end of the surgical session, in order to remove the CENTRAL

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

## 2.4 DEVICE ASSEMBLY


Open the cardboard box containing the TOOL BODY, taking care not to damage the internal part. Do not use scalpels or sharp objects to open the package. Follow all the standard procedures for opening a sterile device. Make certain the device is intact and shows no signs of tampering.


- Open the blister pack and place the contents in a sterile area in the operating field.

- Hold the TOOL BODY container with one hand; ensure that the personnel assigned to this operation are in sterile gowns. Insert the protective cone (INSERTER BOX), included in the package, onto the rear part of the handpiece so as to protect its sterility during introduction of the CENTRAL UNIT.

## 2.2 SETTINGS SELECTION AND SETUP

### Button **A**: **START**

Pressing this button for 5 seconds activates the STAND-BY mode. Pressing button **A** again increases the speed of rotation which is displayed through the blue LEDs located on the grip.

### Button **B**: **CONTINUOUS Mode**


Selecting this mode, the device rotates the tool clockwise. To stop rotation, press button **B** again.

### Button **C**: **ALTERNATING Mode**

Selecting this mode, the device rotates the tool in alternating directions. To stop rotation, press button **C** again.

## 2.3 INSPECTION AND CLEANING

The CENTRAL UNIT can be cleaned using disinfectants normally used in hospitals and healthcare facilities, although care should be taken not to rub too hard as that could damage the surfaces. Soft brushes or sponges should be used, avoiding abrasive material completely in order to prevent damage. The TOOL BODY is supplied in a sterile package, so it must not be resterilized or reused after its first use.

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
## 1.5 PLACE OF USE

*FRII Power* is designed to be used as an operating room device, since it will be used by surgeons to perform arthroscopic procedures. The CHARGING BASE must be placed on a shelf near the mains power socket.

After the CENTRAL UNIT is inserted inside the TOOL BODY container, the device can be held by the surgeon and used for its intended purposes.

## 2. METHOD OF USE

### 2.1 DESCRIPTION OF THE COMMANDS

- A:** START
- B:** CONTINUOUS Mode
- C:** ALTERNATING Mode
- D:** Flow regulator
- E:** Speed display LEDs
- F:** Malfunction indicator LED

- ❑ For no reason whatsoever should the blade or bur be replaced with the motor running.
- ❑ Any damage to the CENTRAL UNIT may result in the leakage of contaminating material; the device must not be used in this case.
- ❑ Strictly follow the instructions regarding the scheduled operating cycle in order to avoid harm to the patient or damage to the device.
- ❑ In Continuous Rotation mode, do not exceed 12 minutes of continued use. For each 12 minutes of continuous use, switch the device off for 5 minutes. An excessive increase in the temperature of the arthroscopic instrument may cause harm or burns to the patient's tissues with which it comes into contact.
- ❑ Always avoid resting the device on the patient.
- ❑ Use exclusively original arthroscopic instruments (request the catalogue from the manufacturer).
- ❑ The TOOL BODY device designed to be used exclusively with the CENTRAL UNIT is a **single-use device supplied in sterile form**.
- ❑ Always avoid leveraging the device as it may become damaged or harm the patient.
- ❑ Applying too much pressure on the instrument may cause it to break and consequently harm the patient.
- ❑ Note that starting the motor causes the instrument to run at the set rpm with rapid acceleration.
- ❑ Grip the device firmly since the rotation of internal components causes the device to vibrate slightly, which could make it slip out of your grip.
- ❑ Use sterile suction tubes only.
- ❑ Before each use, ensure that the suction channel of the TOOL BODY device is clean.
- ❑ Do NOT clean the device using acetone, ether, freon, petroleum derivatives or other solvents.
- ❑ The **FRII Power** device has been tested in reference to Regulations EN 60601-1 and EN 60601-1-2. If the user connects the device to other instruments which have not been previously validated for combined use according to Regulations EN 60601-1 and EN 60601-1-2, the user must ensure that the coupling between the two devices meets the requirements of the aforesaid Regulations. Otherwise the manufacturer disclaims all responsibility.

**The device must be used exclusively as specified herein.**

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- now the CENTRAL UNIT is ready to be inserted into the container of the FRII POWER device.

**Note!**


The CENTRAL UNIT **IS NOT STERILE** and thus should be handled by personnel who are not in sterile gowns, while the TOOL BODY container is in a sterile package and must therefore be handled by personnel in sterile gowns.

While inserting the CENTRAL UNIT inside the TOOL BODY container, maximum care must be taken to ensure that the sterile and non-sterile parts do not come into contact in order to avoid cross contamination.


## 1.4 WARNINGS AND PRECAUTIONS

- No component of the device may be modified for any reason whatsoever.
- If the delivered components show any signs of evident or suspected damage, do not use the components and notify the manufacturer.
- Maintenance must be carried out by manufacturer-authorized personnel only.
- The device described in this manual is designed to be used by medical personnel with extensive knowledge of the surgical techniques requiring the use of motorized instruments.
- Improper use of the device may cause harm to the patient.
- The device must not be used at the same time or near electrosurgical systems, short-wave or microwave therapy equipment.
- The device may not be used in the presence of flammable anesthetic gases or mixtures.
- The surgeon must ensure that the battery of the device is completely charged, as indicated in this user manual, before using the device on the patient.
- It is recommended that protective glasses be worn by medical personnel and assistants present in the area where the device is being used.

- insert the power supply cable into the back of the charging base;

- press the power switch and move it the position "I" ;

- place the CENTRAL UNIT in one of the four special seats on the CHARGING BASE;


- wait until the GREEN LED on the CHARGING BASE indicates that the CENTRAL UNIT is charged;

- **Charge the device for at least 2 hours before use;**
- **Do not charge the device for more than 24 hours;**

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	Manufacturer
	Medical Device compliant with Directive 93/42/EEC as amended (Directive 2007/47)
	Do not dispose of this product as household waste. Arrange for the re-use or recycling of this product as set forth in Directive 2002/96/EC of the European Parliament and European Union Council on the disposal of waste electrical and electronic equipment.

### 1.3 ASSEMBLY AND CONNECTION OF THE COMPONENTS

**Before assembling the system components, check that there are no signs of damage on the components themselves after removing them from the delivery package. If any damage is found or suspected, notify Medicon immediately using the addresses provided at the end of this document.**

The components of the device must be assembled following the steps below:

- place the CHARGING BASE on a level, solid surface at a height of more than 1.50 meters to ensure easy access;

- check that the rear power switch is on position "0";

## 1.2 DESCRIPTION OF THE SYMBOLS USED

The graphic symbols used on the device and in this manual refer to Regulations UNI CEI EN ISO 15223-1, CEI EN 60601-1, the WEEE Directive and Directive 2002/96/ECC.

SYMBOL	DESCRIPTION
	Product code
	Production lot
	Serial number
	Device with type-B applied parts
	Class II EM device
	Caution: read the attached documents before operating the device. Alternatively, the symbol in question may be used as a "Warning" symbol
	See the User Manual
	Single-use Device
 	Sterilized with Ethylene Oxide
	Latex-free
	Manufacturing date
	Expiration date
	Temperature limits

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<b>FR10XXYY</b>		The tool body is composed of the handpiece and tool (Blade or Bur), which actually performs the function of removing the tissue, and a shell that covers and protects the INTERNAL CENTRAL UNIT. A series of buttons allow the user to select the functions. It also has LEDs that indicate the type and levels of the selected function and a suction tube for liquids	The TOOL BODY is designed for the resection of soft tissue or bone tissue in arthroscopic procedures. The functions are controlled through a button panel which displays the active function and its operating level. It also allows the suction of liquids through the assembled tube.	<b>Sterile SINGLE-USE</b>
Where <b>XX=</b> <b>Tool Ø</b>  <b>YY=</b> <b>Tool Model</b>	TOOL BODY			

REF	Name	Description	Function	Supplied Us
<b>FR20700</b>	REMOVAL BOX	This is a plastic-coated box supplied inside the sterile package of the TOOL BODY	The REMOVAL BOX allows the user to remove the central unit from the TOOL BODY after the surgical procedure	<b>Sterile SINGLE-USE</b>
REF	Name	Description	Function	Supplied Us
<b>FR20800</b>	INSERTER BOX	This is a plastic-coated box supplied inside the sterile package of the TOOL BODY	The INSERTER BOX allows the user to insert the central unit into the TOOL BODY before the surgical procedure	<b>Sterile SINGLE-USE</b>

The system is delivered inside a package containing:

- 1 CHARGING BASE complete with power supply cable;
- 1 INTERNAL CENTRAL UNIT;
- 1 TOOL BODY  **SINGLE-USE Component (STERILE package)**

## FOREWORD


The device described below is designed to be used in an operating room by medical personnel with extensive knowledge of the surgical techniques requiring its use.

It must be used exclusively with instruments and accessories produced by Medicon according to the instructions provided in this manual and on the packages.

The device is indicated for arthroscopic procedures on human joints.

## 1. DESCRIPTION

### 1.1 DESCRIPTION OF THE COMPONENTS AND THEIR FUNCTION

REF	Nome	Descrizione	Funzione	Caratteristiche di fornitura
FR20100	<b>KIT COMPOSE D OF CENTRAL UNIT + CHARGING BASE</b>	The CENTRAL UNIT is composed of an integrated unit dedicated to the system's operating functions. The CHARGING BASE is composed of a support designed to hold up to four CENTRAL UNITS.	It provides the tool with rotation and oscillation functions and controls its operating modes. It guarantees an 1 hour runtime. It checks the charge status of the INTERNAL CENTRAL UNIT and charges it if necessary	<b>Non-Sterile and Re-usable</b>
		 		

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## NOTE:

## CONTENTS

FOREWORD	34
<b>1. DESCRIPTION</b>	<b>34</b>
1.1 DESCRIPTION OF THE COMPONENTS AND THEIR FUNCTION	34
1.2 DESCRIPTION OF THE SYMBOLS USED	36
1.3 ASSEMBLY AND CONNECTION OF THE COMPONENTS	37
1.4 WARNINGS AND PRECAUTIONS	39
1.5 PLACE OF USE	41
<b>2. METHOD OF USE</b>	<b>41</b>
2.1 DESCRIPTION OF THE COMMANDS	41
2.2 SETTINGS SELECTION AND SETUP	42
2.3 INSPECTION AND CLEANING	42
2.4 DEVICE ASSEMBLY	44
2.5 OPERATING PARAMETERS	46
2.6 SUCTION	47
2.7 INSTRUCTIONS FOR USE OF THE DEVICE	47
<b>3. CLEANING AND DISINFECTION</b>	<b>49</b>
<b>4. MAINTENANCE</b>	<b>50</b>
<b>5. PRESERVATION AND STORAGE</b>	<b>50</b>
<b>6. TECHNICAL INFORMATION</b>	<b>50</b>
<b>7. SERVICE LIFE OF THE INSTRUMENT</b>	<b>51</b>
<b>8. DISPOSAL</b>	<b>52</b>
<b>9. GUARANTEE</b>	<b>52</b>
<b>10. TROUBLESHOOTING</b>	<b>53</b>
<b>11. INFORMATION PROVIDED BY THE MANUFACTURER REGARDING THE ELECTROMAGNETIC COMPATIBILITY OF THE DEVICE</b>	<b>53</b>

**NOTE:**

Manufacturing by:

MEDICON Medical & Consulting  
Via Moravia 3 - 35030 Selvazzano (PD) - Italy  
Tel. +39 049 8974421 - Fax +39 049 8974421  
e-mail: info@mediconitaly.com

Carefully read the user manual before using FRII Power

**CE** 0476 - Product compliant with Directive **93/42/EEC** as amended (Directive 2007/47) for medical devices, and in compliance with reference standards **CEI EN 60601-1, 60601-1-2, 60601-1-6, CEI IEC 62304 and Directive 2002/96/EEC**

**NOTE:**

## **Motorized system for arthroscopic surgery**

### **User Manual**