

QUADRIL
HIP / CADERA



JOELHO
KNEE / RODILLA



MEDICINA ESPORTIVA
SPORTS MEDICINE / MEDICINA DEPORTIVA



COLUNA
SPINE / COLUMNNA



TRAUMA
TRAUMA / TRAUMA



ICONACY Orthopedic Implants



QUADRIL

HIP / CADERA

**I-TAPPER SYSTEM
I-HIP SYSTEM
I-HIP GT SYSTEM**

COMPONENTE ACETABULAR I-HIP CONSTRAINED

ACETABULAR COMPONENT I-HIP CONSTRAINED
COMPONENTE ACETABULAR I-HIP CONSTRAINED



COMPONENTE ACETABULAR I-HIP

Liga de Titânio Ti-6Al-4V
Wrought Ti-6Al-4V alloy
Aleación de Titânio Ti-6Al-4V
ASTM F136
ISO 5832-3

Revestimento Poroso
Porous Coating
Revestimento Poroso
ASTM F1580

COMPONENTE ACETABULAR 3 FUROS I-HIP ACETABULAR COMPONENT 2 HOLES I-HIP COMPONENTE ACETABULAR 2 AGUJEROS I-HIP	REFERÊNCIA PART NUMBER REFERÊNCIA
OD 48 G2	01-134811
OD 50 G3	01-135011
OD 52 G3	01-135211
OD 54 G3	01-135411

LINER + ANEL DE TRAVAMENTO LINER + LOCKING RING INSERTO + ANILLO DE BLOQUEO	REFERÊNCIA PART NUMBER REFERÊNCIA
G2 ID 28	01-080228
G3 ID 32	01-080332





PARAFUSO ÓSSEO

BONE SCREW
TORNILLO



PARAFUSO ÓSSEO

Liga de Titânio Ti-6Al-4V
Wrought Ti-6Al-4V alloy
Aleación de Titânio Ti-6Al-4V
ASTM F136
ISO 5832-3

PARAFUSO ÓSSEO BONE SCREW TORNILLO ÓSSEO	REFERÊN CIA PART NUMBER REFERÊN CIA
6,5 x 20mm	62-016520
6,5 x 25mm	62-016525
6,5 x 30mm	62-016530
6,5 x 35mm	62-016535
6,5 x 40mm	62-016540
6,5 x 45mm	62-016545
6,5 x 50mm	62-016550
6,5 x 55mm	62-016555
6,5 x 60mm	62-016560