

QUADRIL
HIP / CADERA



JOELHO
KNEE / RODILLA



MEDICINA ESPORTIVA
SPORTS MEDICINE / MEDICINA DEPORTIVA



COLUNA
SPINE / COLUMNNA



TRAUMA
TRAUMA / TRAUMA





QUADRIL

HIP / CADERA

I-HIP REVISION
TM REVISION



SISTEMA DUO MOBILITY

DUO MOBILITY SYSTEM
SISTEMA DUO MOBILITY



Liga Conformada de Co-Cr-Mo
Wrought Co-Cr-Mo alloy
Aleación Conformada de Co-Cr-Mo
ASTM F1537
ISO 5832-12



UHMWPE
ASTM F648
ISO 5834-2

LINER EM UHMWPE



LINER EM UHMWPE - CABEÇA FEMORAL

COMPONENTE ACETABULAR CIM DUO MOBILITY CEMENTED ACETABULAR CUP DUO MOBILITY COPA ACETABULAR CEMENTADA DOBLE MOVILIDAD	REFERÊNCIA PART NUMBER REFERÊNCIA
OD 44	62-634438
OD 48	62-634842
OD 52	62-635246
OD 56	62-635650

LINER EM UHMWPE DUO MOBILITY UHMWPE LINER DUO MOBILITY LINER EN UHMWPE DUO MOBILITY	REFERÊNCIA PART NUMBER REFERÊNCIA
OD 38 P FEMORAL 22	62-663822
OD 42 P FEMORAL 28	62-664228
OD 46 P FEMORAL 28	62-664628
OD 50 P FEMORAL 28	62-665028