

QUADRIL  
*HIP / CADERA*



JOELHO  
*KNEE / RODILLA*



MEDICINA ESPORTIVA  
*SPORTS MEDICINE / MEDICINA DEPORTIVA*



COLUNA  
*SPINE / COLUMNNA*



TRAUMA  
*TRAUMA / TRAUMA*



ICONACY Orthopedic Implants



# QUADRIL

HIP / CADERA

I-HIP REVISION  
TM REVISION

# COMPONENTE ACETABULAR TM (TRABECULAR METAL)

ACETABULAR COMPONENT TM (TRABECULAR METAL)

COMPONENTE ACETABULAR TM (TRABECULAR METAL)



Liga de Titânio Ti-6Al-4V  
Ti-6Al-4V Alloy  
Aleación de Titanio  
ASTM F3001

COMPONENTE ACETABULAR TM

COMPONENTE ACETABULAR 3 FUROS TM (TRABECULAR METAL) ACETABULAR COMPONENT 3 HOLES TM (TRABECULAR METAL) COMPONENTE ACETABULAR 3 AGUJEROS TM (TRABECULAR METAL)		REFERÊN CIA PART NUMBER REFERÊN CIA
Liner Pequeno Small Liner Liner Pequeño	OD 44	62-110144
	OD 46	62-110146
	OD 48	62-110148
Liner Médio Medium Liner Liner Medio	OD 50	62-110250
	OD 52	62-110252
Liner Grande Large Liner Liner Grande	OD 54	62-110354
	OD 56	62-110356
	OD 58	62-110358
	OD 60	62-110360
	OD 62	62-110362
	OD 64	62-110364



# SISTEMA DUO MOBILITY

DUO MOBILITY SYSTEM  
SISTEMA DUO MOBILITY



UHMWPE  
ASTM F648  
ISO 5834-2

**LINER EM UHMWPE**



**LINER EM UHMWPE - CABEÇA FEMORAL**



Liga Conformada de Co-Cr-Mo  
*Wrought Co-Cr-Mo alloy*  
*Aleación Conformada de Co-Cr-Mo*  
ASTM F1537  
ISO 5832-12

**LINER METÁLICO**

LINER EM UHMWPE DUO MOBILITY UHMWPE LINER DUO MOBILITY LINER EN UHMWPE DUO MOBILITY		REFERÊN CIA PART NUMBER REFERÊN CIA
OD 28		62-272822
OD 32	Cabeça Femoral Ø22mm Femoral Head Ø22mm Cabeza Femoral Ø22mm	62-273222
OD 36		62-273622
OD 36	Cabeça Femoral Ø28mm Femoral Head Ø28mm Cabeza Femoral Ø28mm	62-273628
OD 40	Cabeça Femoral Ø22mm Femoral Head Ø22mm Cabeza Femoral Ø22mm	62-274022
	Cabeça Femoral Ø28mm Femoral Head Ø28mm Cabeza Femoral Ø28mm	62-274028
	Cabeça Femoral Ø32mm Femoral Head Ø32mm Cabeza Femoral Ø32mm	62-274032

LINER METÁLICO DUO MOBILITY METALIC LINER DUO MOBILITY LINER METÁLICO DUO MOBILITY		REFERÊN CIA PART NUMBER REFERÊN CIA
ID 28	Extra Pequeno Extra Small Extra Pequeño	62-260028
ID 32	Peque no Small Pequeño	62-260132
ID 36	Médi o Mediu m Medio	62-260236
ID 36	Gran de Large Grande	62-260336
ID 40		62-260340



# PARAFUSO ÓSSEO

BONE SCREW  
TORNILLO



PARAFUSO ÓSSEO

Liga de Titânio Ti-6Al-4V  
Wrought Ti-6Al-4V alloy  
Aleación de Titânio Ti-6Al-4V  
ASTM F136  
ISO 5832-3

PARAFUSO ÓSSEO BONE SCREW TORNILLO ÓSSEO	REFERÊNCI A PART NUMBER REFERÊNCIA
6,5 x 20mm	62-016520
6,5 x 25mm	62-016525
6,5 x 30mm	62-016530
6,5 x 35mm	62-016535
6,5 x 40mm	62-016540
6,5 x 45mm	62-016545
6,5 x 50mm	62-016550
6,5 x 55mm	62-016555
6,5 x 60mm	62-016560