



I-Hip Revision Stem

Your patient specific solution

Distributed by

PROMED
QUIRÚRGICOS



I – HIP REVISION STEM

Tallas de Vástago


TALLAS DEL VASTAGO

El vástago REVISION esta compuesto de 12 tamaños de implantes distales: 6 diámetros de 14mm a 24mm, con tamaños que varían en 2mm y dos longitudes para cada diámetro, 140mm y 200mm.

▼ VASTAGOS DE REVISION

| DIAMETRO | LONGITUD |
|----------|----------|
| 14 mm* | 140 mm |
| 14 mm* | 200mm |
| 16 mm | 140 mm |
| 16 mm | 200 mm |
| 18 mm | 140 mm |
| 18 mm | 200 mm |
| 20 mm | 140 mm |
| 20 mm | 200 mm |
| 22 mm | 140 mm |
| 22 mm | 200 mm |
| 24 mm | 140 mm |
| 24 mm | 200 mm |

* Vástagos 14 no aptos para cuellos lateralizados (prueba según ISO 7206-4)



Adaptabilidad a la anatomía utilizando ángulo de 4° proximal.

I – HIP REVISION STEM

Tallas de Vástagos


▼ CUELLOS DE REVISION

ESTANDAR – 135°:

Tabla de offset con uso de cabeza M (+0 mm)

| ALTURA DE CUELLO | -4° | NEUTRAL STD | +4° |
|------------------|------|-------------|------|
| ■ 50 | 32,1 | 35,6 | 39,1 |
| 60 | 31,4 | 35,6 | 39,8 |
| 70 | 30,7 | 35,6 | 40,5 |
| 80 | 30,0 | 35,6 | 41,1 |
| 90 | 29,3 | 35,6 | 41,8 |
| 100 | 28,6 | 35,6 | 42,5 |
| 110 | 27,9 | 35,6 | 43,2 |

■ A petición del usuario



ACD 135°

LATERALIZADO - 132°:

Tabla de offset con uso de cabeza M (+0 mm)

| ALTURA DE CUELLO | -4° | NEUTRAL STD | +4° |
|------------------|------|-------------|------|
| ■ 50 | 36,6 | 40,1 | 43,6 |
| 60 | 35,9 | 40,1 | 44,3 |
| 70 | 35,2 | 40,1 | 45,0 |
| 80 | 34,5 | 40,1 | 45,7 |
| 90 | 33,8 | 40,1 | 46,4 |
| 100 | 33,2 | 40,1 | 47,1 |
| 110 | 32,5 | 40,1 | 47,8 |

■ A petición del Usuario




ACD 132°

I – HIP REVISION STEM

Set Instrumental

▣ 9038.15.000 Set n.1 para Vástago de Revision Femoral




| Ref. | CODIGO | DESCRIPCION | Cant. |
|------|-------------|------------------|-------|
| A1 | 9038.10.010 | Reamer Dia. 12mm | 1 |
| A1 | 9038.10.015 | Reamer Dia. 13mm | 1 |
| A1 | 9038.10.020 | Reamer Dia. 14mm | 1 |
| A1 | 9038.10.025 | Reamer Dia. 15mm | 1 |
| A1 | 9038.10.030 | Reamer Dia. 16mm | 1 |
| A1 | 9038.10.035 | Reamer Dia. 17mm | 1 |
| A1 | 9038.10.040 | Reamer Dia. 18mm | 1 |
| A1 | 9038.10.045 | Reamer Dia. 19mm | 1 |
| A1 | 9038.10.050 | Reamer Dia. 20mm | 1 |
| A1 | 9038.10.055 | Reamer Dia. 21mm | 1 |
| A1 | 9038.10.060 | Reamer Dia. 22mm | 1 |
| A1 | 9038.10.065 | Reamer Dia. 23mm | 1 |
| | 9038.15.950 | Sterilizable Box | 1 |

I – HIP REVISION STEM


Set Instrumental Adicional

9038.16.000 Instrument Set n.2 for Revision Femoral Stem



| Ref. | CODE | DESCRIPTION | Cant. |
|------|-------------|--|-------|
| A2 | 9095.10.110 | Manual Snap Wrench | 1 |
| B2 | 9038.10.150 | Trial Neck H50 with Screw | 1 |
| B2 | 9038.10.160 | Trial Neck H60 with Screw | 1 |
| B2 | 9038.10.170 | Trial Neck H70 with Screw | 1 |
| B2 | 9038.10.180 | Trial Neck H80 with Screw | 1 |
| B2 | 9038.10.190 | Trial Neck H90 with Screw | 1 |
| B2 | 9038.10.200 | Trial Neck H100 with Screw | 1 |
| B2 | 9038.10.210 | Trial Neck H110 with Screw | 1 |
| C2 | 9038.10.310 | Lateralized Trial Neck H50 with Screw | 1 |
| C2 | 9038.10.320 | Lateralized Trial Neck H60 with Screw | 1 |
| C2 | 9038.10.330 | Lateralized Trial Neck H70 with Screw | 1 |
| C2 | 9038.10.340 | Lateralized Trial Neck H80 with Screw | 1 |
| C2 | 9038.10.350 | Lateralized Trial Neck H90 with Screw | 1 |
| C2 | 9038.10.360 | Lateralized Trial Neck H100 with Screw | 1 |
| C2 | 9038.10.370 | Lateralized Trial Neck H110 with Screw | 1 |
| D2 | 9095.10.511 | Trial Head Low Taper 12/14 Dia. 28mm S | 1 |
| D2 | 9095.10.512 | Trial Head Low Taper 12/14 Dia. 28mm M | 1 |
| D2 | 9095.10.513 | Trial Head Low Taper 12/14 Dia. 28mm L | 1 |
| E2 | 9095.10.521 | Trial Head Low Taper 12/14 Dia. 32mm S | 1 |
| E2 | 9095.10.522 | Trial Head Low Taper 12/14 Dia. 32mm M | 1 |
| E2 | 9095.10.523 | Trial Head Low Taper 12/14 Dia. 32mm L | 1 |
| | 9038.16.950 | Sterilizable Box | 1 |

9038.17.000 Instrument Set n. 3 for Revision Femoral Stem



| Ref. | CODE | DESCRIPTION | Cant. |
|------|-------------|-------------------------|-------|
| A3 | 9038.10.100 | Dynamic Impactor | 1 |
| B3 | 9038.10.110 | Manual Impactor | 1 |
| C3 | 9038.10.115 | Neck Reamer Guide | 1 |
| D3 | 9038.10.120 | Neck Reamer | 1 |
| E3 | 9038.10.230 | Neck Impactor-Extractor | 1 |
| F3 | 9038.10.240 | Neck Stopper | 1 |
| G3 | 9038.10.250 | Trial Neck Extractor | 1 |
| H3 | 9095.10.117 | Allen Wrench | 1 |
| | 9038.17.950 | Sterilizable Box | 1 |