

QUADRIL
HIP / CADERA



JOELHO
KNEE / RODILLA



MEDICINA ESPORTIVA
SPORTS MEDICINE / MEDICINA DEPORTIVA



COLUNA
SPINE / COLUMNNA



TRAUMA
TRAUMA / TRAUMA





QUADRIL

HIP / CADERA

I-HIP REVISION
TM REVISION

PLACA ACETABULAR I-TAPPER REVISION

ACETABULAR PLATE I-TAPPER REVISION

PLACA ACETABULAR I-TAPPER REVISIÓN



Titânio não Ligado GR2
Unalloyed titanium GR2
Titanio no Aleado GR2
ASTM F67
ISO 5832-2

PLACA ACETABULAR I-TAPPER

| PLACA ACETABULAR I-TAPPER REVISION ACETABULAR PLATE I-TAPPER REVISION PLACA ACETABULAR I-TAPPER REVISIÓN | | REFERÊN CIA PART NUMBER REFERÊN CIA |
|---|-------|---|
| | OD 50 | 62-200150 |
| | OD 54 | 62-200154 |
| | OD 58 | 62-200158 |
| | OD 62 | 62-200162 |
| | OD 66 | 62-200166 |
| Long a Large Largo | OD 50 | 62-200250 |
| | OD 54 | 62-200254 |
| | OD 56 | 62-200258 |
| | OD 62 | 62-200262 |
| | OD 66 | 62-200266 |



SISTEMA DUO MOBILITY

DUO MOBILITY SYSTEM
SISTEMA DUO MOBILITY



Liga Conformada de Co-Cr-Mo
Wrought Co-Cr-Mo alloy
Aleación Conformada de Co-Cr-Mo
ASTM F1537
ISO 5832-12



UHMWPE
ASTM F648
ISO 5834-2

LINER EM UHMWPE



LINER EM UHMWPE - CABEÇA FEMORAL

| COMPONENTE ACETABULAR CIM DUO MOBILITY CEMENTED ACETABULAR CUP DUO MOBILITY COPA ACETABULAR CEMENTADA DOBLE MOVILIDAD | REFERÊNCIA PART NUMBER REFERÊNCIA |
|---|---|
| OD 44 | 62-634438 |
| OD 48 | 62-634842 |
| OD 52 | 62-635246 |
| OD 56 | 62-635650 |

| LINER EM UHMWPE DUO MOBILITY UHMWPE LINER DUO MOBILITY LINER EN UHMWPE DUO MOBILITY | REFERÊNCIA PART NUMBER REFERÊNCIA |
|---|---|
| OD 38 P FEMORAL 22 | 62-663822 |
| OD 42 P FEMORAL 28 | 62-664228 |
| OD 46 P FEMORAL 28 | 62-664628 |
| OD 50 P FEMORAL 28 | 62-665028 |